



Iris Valve - Repeatability Test

Brightwater, Seattle, WA. : Brown and Caldwell

- Repeatability test at 90% Open
 - Open to 90% - measure 1,800 scfm
 - Closed to 30% for 30 seconds
 - Open to 90% for 30 seconds
 - Measure 1,791 scfm



Richmond, Indiana – Donohue & Associates

- Repeatability test at 63.2% Open
 - Open to 63.2% - measure 164 scfm
 - Close to 40% for 30 seconds
 - Open to 63.2% for 30 seconds
 - Measure 163.4 scfm



Fishers, Indiana – Bowen Engineering

Repeatability test at 94.4% Open

- Open to 94.4% - measure 982 scfm
- Close to 70% for 30 seconds
- Open to 94.4% for 30 seconds
 - Measure 988 scfm



RESULT: 99.3% REPEATABLE!